

## Ion exchanger IE 2500

Pressure resistant ion exchangers are often the most economical way to demineralise water. The ion exchange resin effectively removes inorganics, i.e. anions and cations, from the water passing through. It is ideal for supplying laboratory glassware washers, autoclaves and sterilizers acc. to DIN 285, air conditioners, cooling systems, ultra-pure water systems as well as for general laboratory purposes. Additionally, it can be used for various industrial applications like feed water for the rinsing of printed circuits and other critical equipment, steam boilers and to produce water equivalent to distilled water for batteries, in conformity with DIN VDE 0510.



## Scope of delivery

Pressure resistant stainless-steel cartridge with handles and lance system Filled with high-quality mixed-bed resin

Туре	IE 2500
Flow rate	500 L/h
Capacity at 300µS/cm*	2500 liters
Pure water quality	0,1 – 20 µS/cm
Material	V4A stainless-steel
Pressure resistant max.	10 bar
Connectors	R ¾"
Max. water temperature	40 °C
Dimensions	237 x 490 mm
Weight	17 kg
Order number	10000200

\* The capacity can vary depending on the impurities in the feed water

Quick connections and additional accessories are available on request. Furthermore, the resin inside the cartridges can be customised depending on the application.